

EPOND Part.No	OEM Part.No.	Description
Nebulizers/Spray chambers/accessories		
Quartz Nebulizers		
EP 2050-M030-C2	47929002	MEINHARD® Type C glass nebulizer, 2ml/min uptake
EP 2050-M030-C3	47929003	MEINHARD® Type C glass nebulizer, 3ml/min uptake 30 PSI
EP 2050-M030-K3	47929005	MEINHARD® Type K glass nebulizer, 3ml/min uptake
EP 2050-M050-C05	47929004	MEINHARD® Type C glass nebulizer, 0.5ml/min uptake
EP 2050-M050-C1	47929001	MEINHARD® Type C glass nebulizer, 1ml/min uptake
EP 2050-M050-C3	47929006	MEINHARD® Type C glass nebulizer, 3ml/min uptake 50 PSI
PFA Nebulizer		
ES-2040-28		MicroFlow PFA-ST self-aspirating nebulizer
Poly-Pro Nebulizer		
EP 2550-4040-57		Poly-Pro-ST self-aspirating nebulizer
PEEK Nebulizer		
EP 3150-0720		Burgener PEEK MiraMist nebulizer
Nebulizer cleaner		
PRL02		Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
Cyclonic Spray Chambers		
EP 2016-0550	31030809	Pyrex cyclonic spray chamber, top socket, 4mm drain
EP 2016-0550ES	47929056	Pyrex cyclonic spray chamber with Easy-Seal nebulizer inlet, PFA cap
EP 2015-015NB	31-031-2779	Cyclonic Spray Chambers, no baffle
EP 2016-0560		Pyrex cyclonic spray chamber, 50ml, top socket, ground glass tapered drain
EP 2016-0570	47929046	Pyrex cyclonic spray chamber with water jacket, top socket, 4mm drain
EP 2016-0600	47929049	Pyrex cyclonic spray chamber, 20ml, top socket, 4mm drain
EP 2016-0620	31031479	Pyrex cyclonic spray chamber, axial, taper to top socket, o-ring inlet, 4mm drain
Spray Chamber		
EP 2015-0140	31026996	Pyrex Scott Type spray chamber, ground-glass joint
Torches/shields		
Quartz Outer Torch Tubes		
ES-1952-9001		Quartz demountable outer torch tube with groundstep, 71mm long
ES-2016-0500	31023722	Quartz demountable outer torch tube, radial, 71mm long
ES-2016-0510		Quartz demountable outer torch tube, radial, 82.5mm long
ES-2016-0700	31031099	Quartz demountable outer torch tube, axial, Panorama 105mm long
Quartz Inner Torch Tubes		
EP 2016-0520	31023723	Quartz demountable inner torch tube
Quartz One-Piece Torch		
EP 2016-0630	49010100	Quartz torch, radial, fixed 1.5mm injector
Torch Assembly		
EP 2016-0400	21-356-110	Fully demountable torch assembly
Injectors and accessories		
Alumina Injector		
EP 2016-0580	31021589	Alumina injector, 2.9mm

EP 2016-0590		Quartz Injector Quartz injector, 2mm
Accessories		
EP 2015-0250	21356040	Sheath Gas Tube Pyrex sheath gas tube
EP 2016-0530	31023724	Centering Ring Centering Ring
EP 2016-0460		O-Ring Kit O-ring kit, demountable torch(Qty 9)
EP 2016-2000	46702525	Viton O-Rings Viton o-rings, pack of 4
EP 2052-0300	31-023-549	PTFE Endcap PTFE adapter(endcap), 6mm nebulizer to 35mm spray chamber