

EPOND Part-No.	OEM Part-No.	Teledyne Leeman Description
Nebulizers/Spray chambers/accessories		
Quartz Nebulizers		
EP 2050-SILQ-50-05		MEINHARD® SiIQ+ Ultra high purity quartz nebulizer
EP 2050-MQ30-A2		MEINHARD® Type A concentric quartz nebulizer, 2 ml/min uptake
EP 2050-MQ30-C1		MEINHARD® Type C concentric quartz nebulizer, 1ml/min
Glass Nebulizers		
EP 2050-M030-C2		MEINHARD® Type C glass nebulizer, 2ml/min uptake
EP 2050-M030-K2		MEINHARD® Type K glass nebulizer, 2ml/min uptake
EP 2050-M030-K3		MEINHARD® Type K glass nebulizer, 3ml/min uptake
PFA Nebulizer		
ES-2040-28		MicroFlow PFA-ST self-aspirating nebulizer
Poly-Pro Nebulizer		
EP 2550-4040-28		Poly-Pro self-aspirating nebulizer
pergo Argon Nebulizer Gas Humidifier		
PRG-01		pergo Argon high solids without nebulizer, includes gas line connector kit
PRG-01-MaPlus		pergo Argon high solids including MEINHARD® Type A Plus high performance nebulizer
Nebulizer Cleaner		
PRL02		Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
FC Gas Quick Connect for MEINHARD® Nebulizers		
EP-FG		FG Gas Quick Connect for MEINHARD® nebulizers
Spray Chambers		
EP 2055-0250	120-3715	Pyrex Scott Type spray chamber, radial
EP 2055-0550	318-00049	Pyrex Scott Type spray chamber, axial
Cyclonic Spray Chambers		
EP 2055-0530	120-00475-2	Cyclonic spray chamber, radial, o-rings for 6 mm nebulizer
EP 2055-0540	120-00475-1	Pyrex cyclonic spray chamber, radial, o-rings for 6mm nebulizer
EP 2055-0580	120-00461-2	Pyrex cyclonic spray chamber for Hildebrand Grid, axial
EP 2055-0590	120-00461-1	Cyclonic spray chamber for Hildebrand Grid, with baffle, axial
Torches/shields		
Quartz One-Piece Torches with Ball Joint		
EP 2055-0620	318-00062	Quartz torch, axial, fixed 2.5mm injector, 187mm length
EP 2055-0640	318-00064	Quartz torch, axial, direct-reading echelle, aqueous, 160mm length
EP 2055-0650	318-00065	Quartz torch, axial, dual view direct-reading echelle, aqueous
EP 2055-0660	120-00496-1	Quartz torch, dual view, DURA
EP 2055-0200	120-00336-1	Quartz torch, radial, fixed 1.5mm injector, tubing and fitting not included
EP 2055-0220	120-00336-2	Quartz torch, radial, fixed 1mm injector, tubing and fitting not included
EP 2055-0230	120-00336-3	Quartz torch, axial, direct-reading echelle, aqueous, fixed 2.5mm injector
Quartz Torches, Demountable Injector		
EP 2055-0690	318-00069	Quartz torch, axial, demountable, direct-reading echelle
EP 2055-0700	318-00070	Quartz torch, axial, demountable, dual view direct-reading echelle
EP 2055-0430	318-0001	Quartz torch, radial, demountable

Accessories

F0 Encapsulated Micro Tubing

EP-F0-25 F0 Encapsulated micro tubing, low uptakes 0.1-0.5ml/min, 0.25mm I.D., 70cm micro capillary

F1 Detachable Capillary

EP-F1-80 F1 Detachable capillary, 0.8mm I.D., 70cm long capillary

F2 Standard PFA Integrated Capillary

EP-F2-50 F2 standard integrated capillary for self-aspiration or pumping, 0.5mm I.D., 70 cm capillary

FitKit 3T

EP-F5-50 FitKit F5 for standard MEINHARD® nebulizers

Tubing and Quick Disconnects

EP 2055-1000 PVC tubing and quick disconnects for Leeman torch

Purge Tubes

EP 2055-A660 318-00066 Purge tube, axial, DV Profile

EP 2055-0670 318-00067 Purge tube, radial, DV Profile

EP 2055-0680 318-00068 Purge tube, radial, Profile