

2014		EPOND Parts list	
EPOND N°	OEM P/No.	Description, accessories for Leeman ICP instruments	CHF
2050-0710N		C-Type standard Concentric Glass Nebulizer, 2 ml/min., <i>with NebLink tube</i>	390
2050-0710NT		'Tyfoon' Hi-performance Concentric Nebulizer, <i>with NebLink sample tube</i>	590
2050-M030-C1	TR30-C1	'Meinhard' Concentric glass nebulizer, Type C, high solids	390
2550-2005		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 400 ul/min uptake version	990
2055-0550		Scott-type Coaxial standard glass spray chamber	470
2052-0300		Adaptor, PTFE 35/6-mm for Scott chamber – glass nebulizer	195
2055-0570		'Althea' Cyclonic glass spray chamber, for standard nebulizers with 6-mm diameter	490
2055-0580		Cyclonic spray chamber, for Hildebrand Grid nebulizer	490
2055-2929		HF-resistant 'Erebus' Cyclonic spray chamber	990
2055-0640	318-00064	Aqueous Axial DRE quartz torch	690
2055-0650	318-00065	Aqueous Dual-View standard DRE 1-piece quartz torch	690
2055-0670	318-00067	Purge Tube, Dual-View DRE	195
2055-0680	318-00068	Radial DRE Purge Tube	270
2055-0690	318-00069	Axial-View Demountable Torch. Aqueous	770
2055-0760	318-00076	Axial-View Torch, Aqueous with injector ID 1,5mm	690
EPOND N°		Flow Meters, Desolvators and Other Accessories	CHF
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply	4970
8100-2051	PRL-02	'PEARL' Cleaning Device for concentric glass nebulizers – <i>suits all types !</i>	215
PERG-001	PRG-01	'PERGO' Argon Humidifier, with tubing and connectors - <i>prevents salting up of nebulizers</i>	990
www. e-pond. biz			