

2015		EPOD Parts List for the Spectro CIROS, GENESIS and ARCOS	
EPOD P/N°	OEM P/No.	Component descriptions - EOP	CHF
EP-2050-0700N	76060510	C-Type Concentric Glass Nebulizer, with 'NebLink' sample tube	390
EP-2050-0700NQ		Standard Concentric Quartz Nebulizer, with 'NebLink' sample tube	660
EP-2050-0710NT		High-Performance 'Tyfoon' Concentric Glass Nebulizer, with sample tube	570
2050-M030-C1	TR-30-C1	'Meinhard' Concentric glass nebulizer, Type C	390
EP-2550-2002	ES-2002-0002	HF-resistant PFA MicroFlow Concentric Nebulizer, 100 ul/min uptake version	990
EP-2550-2005	ES-2005-0002	HF-resistant PFA MicroFlow Concentric Nebulizer, 400 ul/min uptake version	990
EP-2550-2028	ME4040	PolyPro-ST HF-res concentric nebulizer, with detachable sample uptake tube	590
EP-2550-2024	ME2040	PFA-ST HF-res concentric nebulizer, with detachable sample uptake tube	990
EP-2030-2800	ES-2030-28	PFA-HS High-Solids HF-res concentric nebulizer, with detachable sample uptake tube	990
EP-2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
EP-2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
EP-5017-0750-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,8 mm	290
EP-5017-0500-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,5 mm	290
EP-3150-0710	TMM3500	Teflon Burgener MiraMist Parallel-Path Nebulizer	990
EP-3150-0720	PMM4000	PEEK Burgener MiraMist Parallel-Path Nebulizer	990
EP-2050-092S	75060502	Fixed Cross-Flow for Spectro, 34mm Inner Diameter	970
EP-2051-020S	48205036	Spectro Modified Lichte Nebulizer	990
EP-2075-0120	48105003	Spectro Small Conical Spray Chamber, with Impact Sphere	550
EP-2075-0130	48105051	Spectro Glass Spray Chamber, Double Pass 90° Exit	550
EP-2075-0130	4810501	Spectro Quartz Spray Chamber, Double Pass 90° Exit	890
EP-3167-1001	ES-3167-1001	Quartz cyclonic spray chamber, can be used with PC3 Peltier cooling unit	890
EP-2075-0140	48105006	Spectro Large Double Pass Spray Chamber, RS12/5 outlet	550
EP-2100-2477		HF-resistant PFA cyclonic spray chamber	970
EP-2256-0470		47-mm diameter PFA spray chamber for ICP-OES, with 'PureCap' end-cap	1590
EP-2256-5470		47-mm PFA spray chamber for Spectro ICP-MS, for use with ESI torch	1190
EP-2100-2478		PolyPropylene adaptor for use of Spectro X-Flow nebulizer with 47-mm sp. Chamber	255
EP-4157-1000		PC3 Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply	2890
EP-9100-5700		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3980
EP-2250-0470	ES-2250-0470	PFA Sapphire HF-resistant Kit for Spectro CIROS, with PFA-ST nebulizer	3990
EP-2075-0170	48105062	Spectro Elbow Adapter, EOP torch to Chamber	70
EP-2075-0180	48105002	Spectro Elbow Adapter Slurry, EOP Torch to Chamber	70
EP-2075-0300		Spectro Drain/Pressure Lock	215
EP-2075-0350		Spectro U-Drain Trap, 13mm OD	70
EP-2075-1702	ES-1702-0000	Quartz Torch for demountable <i>axial-view</i> CIROS	390
EP-2075-0210	48105008	Spectro Quartz Injector/Aerosol Tube 1.5 mm	155
EP-2075-0220	48105010	Spectro Quartz Injector/Aerosol Tube 1.8 mm	155
EP-2075-0230	48105009	Spectro Quartz Injector/Aerosol Tube 2.5	155
EP-1980-0180	ES-1980-0180	<i>Sapphire Injector in PFA adaptor, ID = 1,8 mm</i>	890
EP-1980-0300	ES-1980-0300	<i>Sapphire Injector in PFA adaptor, ID = 3 mm</i>	890
EP-1981-0200	ES-1981-0200	Platinum Injector in PFA adaptor, ID = 2 mm	1790
EP-2075-0400	48105052	Spectro EOP Flared End Quartz Torch 2.5mm injector, 4mm side arms	390
EP-2075-0400F	75160526	Spectro EOP Flared End Quartz Torch 2.5mm injector, with fittings	390
EP-2075-0400LL		Spectro EOP Flared End Long-Life Quartz Torch	490
EP-2075-0410	48105052A	Spectro EOP Fixed Flared End Quartz Torch, 1.8mm ID injector	390
EP-2075-0420	48206007	Spectro EOP Demountable Torch Tube Assembly, ML-175004	330
EP-2075-0430	75160538	PTFE Housing for Demountable Torch	270
EP-2075-0440		Spectro EOP Demountable Flared End Torch	490
EP-2075-0450	75160545	Complete Demountable Torch	850
EP-2075-0460	48105058	Quartz Injector Tube, 2.5mm ID	190

EP-2075-0480	75160537	Spectro EOP Fixed Sheath Gas System Torch, OEM 75160549	550
EP-2075-0550		Gas Connectors with Part Nos. 44302072 and 443302077	50
EP-2075-0930		Spectro Replacement O-Ring for Endcap	8
EP-1902-0210		1-piece Quartz Torch for radial-view CIROS, with 1,0mm injector	500
EP-2075-1903		1-piece Quartz Torch for radial-view CIROS, with 1,8mm injector	500
EP-1902-0225		1-piece Quartz Torch for radial-view CIROS, with 2,5mm injector	500
EP-2075-1992		Complete Demountable Torch with quartz injector, for <i>radial-view</i> Spectro CIROS	810
EP-2075-1902	ES-1902-0000	Quartz Torch for demountable <i>radial-view</i> CIROS torch P/No. EP-2075-1992	390
EP-1983-2102		Quartz Injector with 1,0mm ID, for use with EP-2075-1992 Quartz radial Torch	420
EP-1983-2152		Quartz Injector with 1,5mm ID, for use with EP-2075-1992 Quartz radial Torch	420
EP-1983-2202		Quartz Injector with 2,0mm ID, for use with EP-2075-1992 Quartz radial Torch	420
EP-1980-00SA		Sapphire Injector for use with EP-2075-1992 Quartz Torch, ID range 1,0 – 2,0 mm.	890
EP-1981-00PT		Platinum Injector for use with EP-2075-1992 Quartz Torch, ID range 1,0 – 2,0 mm.	2150
EPOND No.		Desolvators and other accessories	CHF
APEX-E-SP		APEX-E Desolvator, with PolyPro-ST nebulizer and ' PLUTO ' Flow Meter	POA
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply	4890
EP-8100-2051	PRL_03	'PEARL' Cleaning Device for concentric glass nebulizers - for all types !	195
EP-2045-0700	PRG-01	'PERGO' Argon Humidifier, with tubing and connectors	990
EP-5016-3250	ES-5016-3250	Spectro AS400 Autosampler Probe, 0,25-mm ID, ~ 100 ul/min.	300
EPOND No.	OEM P/No.	Descriptions of tubing accessories	CHF
<i>See web site</i>		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12	48
<i>See web site</i>		<i>Cone-terminated</i> Peri-Pump tubes, for low-uptake apps., pkt of 10	63
ERT-040		EPOND 'RubeTubes' , Tygon-core sample tube connectors, packet of 40	165
DLK-075		<i>DrainLink spray chamber drain</i> adaptor and tube with 0,75mm ID, pkt 5	52
		Autosamplers and ICP / AAS Standard and Reference Solutions	
SC-8000-XXXX		<i>Complete range of ESI SC-series autosamplers - please ask for details</i>	POA
ICP-standards		<i>Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards</i>	
		... custom standards supplied - according to your specifications ...	
		www.e-pond.biz	