

| 2015 | | SpectroFlame, Spectro Modula and SOP instruments | |
|--------------|-----------|---|------|
| EPOND P/N° | OEM P/No. | Description | CHF |
| 2050-0710N | | K-Type Concentric Glass Nebulizer, with 'NebLink' sample tube | 390 |
| 2050-M030-C1 | TR-30-C1 | ' Meinhard ' Concentric glass nebulizer, Type C | 390 |
| 2050-0710NT | | High-Performance ' Tyfoon ' Concentric Glass Nebulizer, with 'NebLink' sample tube | 570 |
| 2550-2005 | | HF-resistant MicroFlow PFA Concentric Nebulizer, 360 ul/min uptake version | 990 |
| 2550-2028 | ME4040 | PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube | 590 |
| 2550-2024 | ME2040 | PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube | 990 |
| 2550-2042 | ES-2042 | Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min. | 90 |
| 2550-2041 | ES-2041 | Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min. | 90 |
| 3150-0710 | TMM3500 | Teflon MiraMist Parallel-Path Nebulizer | 990 |
| 3150-0720 | PMM4000 | PEEK MiraMist Parallel-Path Nebulizer | 990 |
| 2051-0100N | | V-Groove High-Solids Glass Nebulizer, with <i>NebLink</i> | 780 |
| 3150-0730 | | Teflon T2100 'Work-horse' Parallel-Path Nebulizer | 990 |
| 2050-092S | | Fixed Cross-Flow for Spectro, 34mm Inner Diameter | 970 |
| 2051-020S | | Spectro Modified Lichte Nebulizer | 990 |
| | | | |
| 2075-0120 | 48105003 | Spectro Small Conical Spray Chamber, with Impact Sphere | 490 |
| 2075-0130 | 48105051 | Spectro Glass Spray Chamber, Double Pass 90° Exit | 490 |
| 2075-0130Q | 48105051Q | Spectro Quartz Spray Chamber, Double Pass 90° Exit | 870 |
| 2075-0140 | 48105006 | Spectro Large Double Pass Spray Chamber, RS12/5 outlet | 490 |
| 2075-0150 | | Spectro Cyclonic Spray Chamber, 12/5 Socket, 13mm Drain | 490 |
| EP-3185-2008 | | Polypropylene Cyclonic spray chamber with RS12/5 socket | 1090 |
| EP-2127-0341 | | PFA adaptor for fitting to RS12/5 spherical joint | 110 |
| EP-5509-0085 | | Custom HF-resistant adaptor, 7-mm injector to RS12/5 | 185 |
| 2075-0160 | | Spectro MLN Cyclonic Spray Chamber | 490 |
| EP-2100-2477 | | HF-resistant PFA spray chamber kit with PP adaptor for Spectro X-Flow nebulizer | 1295 |
| EP-2256-0470 | | 47-mm diameter PFA spray chamber for ICP-OES, with 'PureCap' end-cap | 1785 |
| EP-2256-5470 | | 47-mm PFA spray chamber for Spectro ICP-MS, for use with ESI torch | 1180 |
| EP-2100-2478 | | PolyPropylene adaptor for use of Spectro X-Flow nebulizer with 47-mm sp. Chamber | 210 |
| EP-4157-1000 | | PC3 Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply | 2870 |
| 2075-0170 | 48105062 | Spectro Elbow Adapter, EOP torch to Chamber | 70 |
| 2075-0180 | 48105002 | Spectro Elbow Adapter Slurry, EOP Torch to Chamber | 70 |
| | | | |
| 2075-0050 | 75060596 | Spectroflame Quartz Torch, 1.8mm injector, 4mm side arms | 490 |
| 2075-005A | 75060596 | Spectro Quartz Torch, 1.8 Injector, with connectors | 490 |
| 2075-005W | | Spectro Wear Metals Quartz Torch, 1.8 Injector | 620 |
| 2075-0060 | 48206002 | Spectro Radial Demountable Torch Tube Assembly | 490 |
| 2075-0070 | 48205001 | Spectro Quartz Torch, 1.0 Injector | 500 |
| 2075-0080 | | Spectro Quartz Torch, 2.5 Injector | 500 |
| 2075-008N | | Spectro (Nitrogen) Quartz Torch | 740 |
| 2075-0009N | | Spectro Demountable Nitrogen Torch | 490 |
| 2075-0210 | 48105008 | Spectro Qtz Injector/Aerosol Tube 1.5 mm | 170 |
| 2075-0220 | 48105010 | Spectro Qtz Injector/Aerosol Tube 1.8 mm | 170 |
| 2075-0230 | 48105009 | Spectro Qtz Injector/Aerosol Tube 2.5 | 170 |
| 2075-0370 | | Spectro Quartz 12/5 Socket, L=35mm | 30 |
| 2075-0380 | | Spectro Quartz 12/5 Ball, L=35mm | 30 |
| 2075-0430 | 75160538 | PTFE Housing for Demountable Torch | 270 |
| 2075-0450 | 75160545 | Complete Demountable Torch | 770 |
| 2075-0460 | 48105058 | Quartz Injector Tube, 2.5mm ID | 190 |
| 2075-0500 | | Spectro Replacement O-Ring | 8 |

| | | | |
|-----------|--|-------------------------------------|----|
| 2075-0550 | | Gas Connectors 44302072 & 443302077 | 50 |
|-----------|--|-------------------------------------|----|

| EPOND No. | | Accessory Description | CHF |
|--------------|--------------|---|------|
| EP-4599-3002 | ES-4599-300A | Spiro TMD Membrane Desolvation Module (latest model of ACM), with Power Supply | 3990 |
| 8100-2051 | PRL-02 | 'PEARL' Cleaning Device for concentric glass nebulizers | 195 |
| 2045-0700 | PRG-01 | 'PERGO' Argon Humidifier, new technology, uses vapour permeation | 990 |
| EP-5016-3150 | ES-5016-3150 | Spectro AS400 Autosampler Probe, 0.15-mm ID, ~ 20 ul/min. | 335 |
| EP-5016-3200 | ES-5016-3200 | Spectro AS400 Autosampler Probe, 0.20-mm ID, ~ 50 ul/min. | 335 |
| EP-5016-3250 | ES-5016-3250 | Spectro AS400 Autosampler Probe, 0.25-mm ID, ~ 100 ul/min. | 335 |

| EPOND No. | OEM P/No. | Description of tubing accessory | CHF |
|---------------------|-----------|---|-----|
| <i>See web site</i> | | All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12 | 48 |
| ERT-040 | | EPOND 'RubeTubes', sample tube connectors, packet of 40 | 165 |
| EST-050 | | PTFE Sample Tube, OD=1,3mm ID=0,50mm L=10m. | 65 |
| DLK-075 | | <i>DrainLink spray chamber drain</i> adaptor and tube with 0,75mm ID, pkt 5 | 52 |

| | | | |
|--|--|--|--|
| Autosamplers and ICP / AAS Standard and Reference Solutions | | | |
|--|--|--|--|

| | | | |
|---------------|--|---|-----|
| SC-8000-XXXX | | <i>Complete range of ESI SC-series autosamplers - please ask for details</i> | POA |
| ICP-standards | | <i>Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards</i> | |

| | | | |
|---|--|--|--|
| www.e-pond.biz | | | |
|---|--|--|--|