

2015	EPOND Parts List : Thermo iCAP Q ICP-MS		
EPOND N°	OEM P/No.	Description	CHF
TRP-50-A0.5	1161870	Low-uptake concentric glass nebulizer, 200uL/min.	390
8100-2051	PRL-02	'Pearl' cleaning device for glass and quartz concentric nebulizers	190
ME2050-52	1162030	PFA-LC Concentric PFA Nebulizer, with gas and sample tubes	990
ES-2040-73	1317090	PFA-ST Concentric PFA nebulizer, ME2040-73	990
EME2020	1600427	Concentric PFA MicroFlow nebulizer, 20 uL/min, ES-2020-0002	990
EME2000	1600165	Concentric PFA MicroFlow nebulizer, 50 uL/min, ES-2000-0002	990
EME2002	1600342	Concentric PFA MicroFlow nebulizer, 100 uL/min, ES-2002-0002	990
2550-2028	ES-4040-28	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube, ME4040	570
2550-2024	ES-2040-27	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube, ME2040	990
ES-2045		Detachable sample capillary tube for PFA-ST nebulizer, ID = 0,15mm., nominal uptake 20 ul/min.	90
ES-2043		Detachable sample capillary tube for PFA-ST nebulizer, ID = 0,20mm., nominal uptake 50 ul/min.	90
ES-2042	1335460	Detachable sample capillary tube for PFA-ST nebulizer, ID = 0,25mm., uptake 100 ul/min.	90
ES-2041	1335470	Detachable sample capillary tube for PFA-ST nebulizer, ID = 0,50mm., uptake 400 ul/min.	90
ES-2050-52	1162030	PFA-LC Concentric PFA nebulizer, includes adapter for LC-ICP-MS and IC-ICP-MS, ME2050-52	990
ME2000	1600165	PFA MicroFlow nebulizer, 50 uL/min	990
ME2002	1600342	PFA MicroFlow nebulizer, 100 uL/min	990
ME2020	1600427	PFA MicroFlow nebulizer for very low uptakes, 20 uL/min	990
PMM4000	4600356	Burgener MiraMist, High-Solids nebulizer (0,4 to 2 mL/min)	990
EAM5000	4600495	Burgener AriMist, Low-uptake, High-Solids, HF-resistant nebulizer (0,2 to 2 mL/min)	990
2046-0001		Obstruction Removal Kit for ST-type nebulizers, with fused silica capillaries	95
Spray Chambers			
ES-3173-1111-23	1317080	Quartz cyclonic spray chamber, ML148030	595
ES-3173-3011-21	1320260	PFA cyclonic spray chamber, ME148032	1190
ES-3150-0073-L	1320270	PFA elbow with add gas port, ME31500073L	155
ES-4190-60-29		PC3 Sample Introduction Kit for iCAP-Q with Quartz Spray Chamber and PFA-ST nebulizer	3385
ES-9190-60-29		PC3-LT Sample Introduction Kit for iCAP-Q with Quartz Spray Chamber and PFA-ST nebulizer	4580
ES-2502-1000	1305780	Additional Gas Line	80
ML148010	1230790	iCAP Q Quartz Torch, ML148010	190
MC31000	1231000	iCAP Q torch holder for quartz torch 1230790	690
ME148017	1305650	1.5mm quartz injector	220
ME148016	1305640	2.0mm quartz injector	220
ME148017	1305600	2.5mm quartz injector	220
ME148050	1323600	1.8mm ID sapphire injector, I73-S18	725
ME148018	1321200	2.0mm platinum injector	1595
OEM P/No. Cones, Multiplier and RF Load Coil			
ES-3000-1221-C	3600812	Nickel Sampler Cone	300
MC11870	1311870	Ni Skimmer Cone, insert ready	455
ES-3000-1226	3601289	Pt Sampler Cone, MC01289	1690
MC24540	1324540	Pt Skimmer Cone, insert ready	1890
MC11880	1311880	Skimmer cone, insert 2.8, High sensitivity	155
MC18470	1318470	Skimmer cone insert 3.0, Standard for A and C, High matrix	155
MC18480	1318480	Skimmer cone insert 3.5, High Matrix Tolerance	155
ES-3000-1088_10	1310900	Sampling cone gasket, pkt 10, MC10900	30
EMC9-9710	1299710	Electron Multiplier (SEM) iCAP Q	2340
EMC0-5260	1305260	RF Load Coil assembly iCAP Q	345
EPOND No.	OEM P/No.	Autosamplers, components and probes	CHF

SC-2, SC-4, etc.		'SC-series' autosamplers are available – to suit any applications – please ask for details	POA
ES-5016-3200-080	ES-5016-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, PEEK support, ~ 50 ul/min., L = 80cm	245
ES-5016-3250-100	ES-5016-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, PEEK support, ~ 100 ul/min., L = 100cm	245
ES-5016-3500-120	ES-5016-3500-120	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, PEEK support, ~ 400 ul/min., L = 120cm	245
Rotor-PF-P6R	PF-P6R	PFA Rotor (1mm bore) for P6 and P7 FAST valves	82
ES-5037-3505-150	1317120	0.5mm ID Carbon fiber SC AS Probe	225
ES-5036-3200-080	ES-5036-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L = 80cm	225
ES-5036-3250-100	ES-5036-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L = 100cm	225
ES-5036-3500-120	ES-5036-3500-120	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, Carbon Fibre support, ~ 400 ul/min., L = 120cm	225
ES-5027-3200-100	ES-5027-3200-100	SC-autosampler probe without flange , ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L=100cm	225
ES-5027-3250-100	ES-5027-3250-100	SC-autosampler probe without flange , ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L=100cm	225
ES-5027-3500-100	ES-5027-3500-100	SC-autosampler probe without flange , ID=0.50mm, Carbon Fibre support, ~ 400 ul/min., L=100cm	225
ES-5056-3200-080	ES-5056-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Ultem support, ~ 50 ul/min., L = 80cm	245
ES-5056-3250-100	ES-5056-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Ultem support, ~ 100 ul/min., L = 100cm	245
ES-5056-3500-100	ES-5056-3500-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, Ultem support, ~ 400 ul/min., L = 100cm	245
ES-5037-3205-100	ES-5037-3205-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L=100cm	225
ES-5037-3255-100	ES-5037-3255-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L=100cm	225
ES-5037-3505-100	ES-5037-3505-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.5mm, Carbon Fibre support, ~ 400 ul/min., L=100cm	225
ES-5037-3505-150	1317120	0.5mm ID Carbon fiber SC AS Probe	225
SC-5037-3991-100		FAST sample probe, ID=1mm, capillary L=100cm	225
SC-5037-3995-150		FAST sample probe, ID=1mm, capillary L=150cm	225
EPOND No.	OEM P/No.	Other Accessories	CHF
EP-PERGO	PRG-01	'PERGO' argon humidifier - 10% discount if purchased with a nebulizer.	1099
ES-4599-3002		Spiro TMD Membrane Desolvation Module (latest vsn of ACM), with Power Supply	3990
SC-0400	SC-0400	oneFAST Basic System with FAST valve for non-ESI autosampler, for X-Series2	6990
8100-2051	PRL-02	'PEARL' Cleaning Device for all types of concentric glass nebulizers	190
EPOND N°	OEM P/No.	Description, tubes	CHF
PVCQ-10-3F	1329400	PVC Flared 3-stop peri-pump tubing, orange/blue, ID=0,25mm, pack of 12	50
PVCQ-15-3F	1320040	PVC Flared 3-stop peri-pump tubing, orange/green, ID=0,38mm, pack of 12	50
PVCQ-20-3F	1320050	PVC Flared 3-stop peri-pump tubing, orange/yellow, ID=0,5mm, pack of 12	50
PVCQ-30-3F	1320060	PVC Flared 3-stop peri-pump tubing, black/black, ID=0,76mm, pack of 12	50
SANQ-51-3	1320110	Santoprene 3-stop peri-pump tubing, Grey/Grey, ID=1,3mm, pack of 12	95
SANQ-125-3	1320120	Santoprene 3-stop peri-pump tubing, Black/White, ID=3,18mm, pack of 12	95
SOLQ-10-3F		Solvent-Flex Tubing, Orange/Black, 3-stops, flared ends, ID=0,25mm, pack of 12	62
SOLQ-15-35		Solvent-Flex Tubing, Orange/Green, 3-stops, flared ends, ID=0,38mm, pack of 12	62
		Autosamplers and ICP / AAS Standard and Reference Solutions	CHF
SC-8000-XXXX		Complete range of ESI SC-series autosamplers and spare parts - please ask for details	POA
TFE-042	1041071	Teflon capillary, OD/ID = 1,6/1,0mm, L=3 metres	25
ME2502-1000	1305780	Addition gas line	40
ME2502-20-R2	1317100	20cm PFA argon line	80
MP500QD	1320280	Spray chamber drain line	15
ME2042	1335460	Nebulizer sample capillary, 100uL/min	80
ME2041	1335470	Nebulizer sample capillary, 400uL/min	80
MM1548	1600061	Teflon capillary, OD/ID = 1,6/0,5mm, L=3 metres	35
MP2-4CH		Peristaltic Pump, 4-channels, autosampler peripheral	2390
ICP-standards		... custom standards supplied - according to your specifications ...	
		www.e-pond.biz	