

2015	EPOND Accessories List - Thermo-Fisher ICP-OES instruments		
Epond N°	OEM P/No.	ICAP-6300, 6500, IRIS Intrepid, Advantage, Axial + DUO	CHF
2050-0800N	13436500	Concentric Glass Nebulizer, K-type, <i>with NebLink sample tube</i>	390
TQ-30-A3		Concentric Quartz Nebulizer, Type A	590
2050-0710N		Standard Concentric Glass Nebulizer, for iCAP series	390
2050-0710NT	14186300	Tyfoon High-performance concentric glass nebulizer, with NebLink capillary	530
TR-30-K2	8423 120 51431	Standard Concentric Glass Nebulizer, for iCAP series, 13960000	390
TR-50-C2	8423 120 51331	<i>Meinhard Type-C</i> Concentric Glass Nebulizer for iCAP series, high solids	530
2051-0100QD	8423 120 51321	V-Groove nebulizer, for iCAP series	870
2550-2028	ME4040	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	590
2550-2024	ME2040	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	990
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
3150-0730	13638800	Burgener PTFE T2002 'High-Solids' Parallel-Path Nebulizer	990
TMM3500	8423 120 51471	Teflon MiraMist Burgener Parallel-Path Nebulizer, for iCAP series	990
BRC-203X		Polyethylene capillaries for Burgener nebulizer, pkt 10	90
BRG-204		Upchurch gas adaptor and argon tube for Burgener nebulizer	30
BRC-20		Upchurch gas adaptor for Burgener nebulizer	10
2052-0300		Teflon Spray Chamber End Cap, for Open-End Spray Chambers, 6mm nebs	290
2027-5150	13462000	TJA Axial / DUO Trace Cyclonic Glass Spray Chamber, 12mm drain	490
2027-515NB		TJA Axial / DUO, Glass Spray Chamber like 2027-5150, No Baffle	490
2027-5160	13612000	TJA Axial / DUO ' Althea ' Cyclonic Glass Spray Chamber, with elbow	490
2027-5163		TJA Axial / DUO ' Libra ' Baffled Cyclonic Spray Chamber	510
2027-5166		TJA Axial / DUO ' Electron ' Mini-Cyclonic Glass Spray Chamber	490
2027-5169		TJA Axial / DUO ' Positron ' Baffled Mini-Cyclonic Spray Chamber	510
2028-5150	8423 120 51411	Glass Cyclonic Spray Chamber, iCAP 6000 series	270
2028-5130	8423 120 51311	Glass Cyclonic Spray Chamber with baffle, iCAP 6000 series	270
2028-5135	8423 120 51441	Jacketed Glass Cyclonic Spray Chamber with baffle, iCAP 6000 series	490
2028-5260	8423 120 51251	Spray Chamber Adaptor for iCAP 6000 DUO series	110
2027-2929		EPOND ' Erebus ' HF-resistant Cyclonic spray chamber	995
EP-4168-1000	ES-4168-1000	PC3 Peltier-cooled Cyclonic Quartz Spray Chamber, with nebulizer and Power Supply	2890
EP-9100-6800		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3990
2025-0230	13768900	TJA DUO Quartz Sleeve with Optical Port, ID = 22,5 mm minimum	165
2027-5100	126432-03	TJA Axial, Trace Torch with Quartz Aqueous Injector, ID = 1,5mm	190
2027-510SI	126440-03	TJA Axial, Trace Torch without Injector, (EP-1504-0200)	570
2027-510DV	138562-00	TJA DUO, Trace DV Torch with 1.5 Aqueous Injector	830
2027-510DVSI	126440-05	TJA DUO, Trace DV Torch without Injector	660
2027-510HF	138566-03	TJA Axial, Trace Torch 1.5 Alumina Injector	1150
EP-1423-0100		Quartz 'Drop-In' Injector with ID = 1,0 mm	190
EP-1423-0150		Quartz 'Drop-In' Injector with ID = 1,5 mm	190
EP-1403-0100		Sapphire 'Drop-In' Injector with ID = 1,0 mm	770
EP-1403-0180		Sapphire 'Drop-In' Injector with ID = 1,8 mm	770
EP-1413-0180		Platinum 'Drop-In' Injector with ID = 1,8 mm	1590
EP-1413-0200		Platinum 'Drop-In' Injector with ID = 2,0 mm	1590
EP-2230-1470	ES-2230-1470-00	ESI Torch with sapphire injector	995
2025-0400	125407-01	TJA Radial, Axial and DV Quartz 1.0 Injector for Organics	235
2025-0410	125407-00	TJA Radial, Axial and DV Quartz 1.5 Injector for Aqueous	235
2025-0420	136502-01	TJA Radial, Axial and DV Quartz 2.0 Injector	235
2025-0430	135602-02	TJA Radial, Axial and DV Quartz 2.4 Injector	235
2027-5250	085207-00	TJA Torch Anti Arc Cap	85
2027-5260	085207-01	TJA Torch Anti Arc Cap / Flat Edge	85
2028-5100	8423 120 51241	iCAP 6000 series DUO Quartz Torch, older type	270
2028-5400	8423 120 51341	iCAP 6000 series Centre Tube, 1.0mm	130
2028-5410	8423 120 51351	iCAP 6000 series Centre Tube, 1.5mm	130
2028-5420	8423 120 51371	iCAP 6000 series Centre Tube, 2.0mm, L = 97mm	130
2028-5270	8423 120 51591	Duo Periscope Quartz Window/Sleeve – ML-128-073	150

2028-5120	8423 120 51841	iCAP 6000 series DUO 'EMT' Quartz Torch, ML128050	370
2028-8072	8423 120 51971	iCAP 6000 series Centre Tube for EMT torch, 2.0mm, ML128072, L = 104,2mm	140
EPOND No.		Autosamplers and probes	CHF
SC-2, SC-4, etc.		'SC-series' autosamplers are available – to suit any applications – please ask for details	POA
ES-5016-3200-080	ES-5016-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, PEEK support, ~ 50 ul/min., L = 80cm	245
ES-5016-3250-100	ES-5016-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, PEEK support, ~ 100 ul/min., L = 100cm	245
ES-5016-3500-120	ES-5016-3500-120	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, PEEK support, ~ 400 ul/min., L = 120cm	245
ES-5036-3200-080	ES-5036-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L = 80cm	245
ES-5036-3250-100	ES-5036-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L = 100cm	245
ES-5036-3500-120	ES-5036-3500-120	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, Carbon Fibre support, ~ 400 ul/min., L = 120cm	245
ES-5027-3200-100	ES-5027-3200-100	SC-autosampler probe without flange , ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L=100cm	225
ES-5027-3250-100	ES-5027-3250-100	SC-autosampler probe without flange , ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L=100cm	225
ES-5027-3500-100	ES-5027-3500-100	SC-autosampler probe without flange , ID=0.50mm, Carbon Fibre support, ~ 400 ul/min., L=100cm	225
ES-5056-3200-080	ES-5056-3200-080	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Ultem support, ~ 50 ul/min., L = 80cm	245
ES-5056-3250-100	ES-5056-3250-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Ultem support, ~ 100 ul/min., L = 100cm	245
ES-5056-3500-100	ES-5056-3500-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.50mm, Ultem support, ~ 400 ul/min., L = 100cm	245
ES-5037-3205-100	ES-5037-3205-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.20mm, Carbon Fibre support, ~ 50 ul/min., L=100cm	245
ES-5037-3255-100	ES-5037-3255-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.25mm, Carbon Fibre support, ~ 100 ul/min., L=100cm	245
ES-5037-3505-100	ES-5037-3505-100	SC-autosampler probe flanged for -ST nebulizer, ID=0.5mm, Carbon Fibre support, ~ 400 ul/min., L=100cm	245
EPOND No.		Accessory Description	CHF
EP-PERGO	PRG-01	'PERGO' argon humidifier – 10% discount if purchased with a nebulizer.	1090
SC-5037-3991-100	"	FAST sample probe, ID=1mm, capillary L=100cm	245
SC-5037-3995-150		FAST sample probe, ID=1mm, capillary L=150cm	245
8100-2051	PRL-02	'PEARL' Cleaning Device for ALL TYPES of concentric glass nebulizers	215
EPOND Part No	OEM P/No.	Description, pumps and tubing	CHF
EPT0007	8423-120-51511, x 2	Aqueous sample Tygon peri-pump tubes, ID=0,64mm., packet of 12, equiv. ES-4398-1064, org-white	40
EPT0012	8423-120-51521, x 2	Aqueous drain Tygon peri-pump tubes, ID=1,02mm., packet of 12, equiv. ES-4398-1102, white-white	40
See web site		'CT – style' Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 12	40
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40	165
SC-8000-XXXX		Complete range of ESI SC-series autosamplers and spare parts - please ask for details	POA
SC-0305-1		Autosampler 5-litre rinse bottle, for SC-DX autosamplers	155
LR-60-16		Rack, 60 (5x12) positions, for 16mm OD vials	40
MP2-4CH		High-stability Peristaltic Pump , 4-channels, autosampler peripheral	2370
EPME-PG41-100		100-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	390
EPME-PG41-500		500-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	790
		reference solutions made to measure - on request	
www. e-pond. biz			