

2015	EPOND Switzerland ICP Parts List		
EPOND N°	OEM P/No.	Thermo IRIS Advantage Radial, Intrepid Radial, and ICAP radial-view	CHF
2050-0800N	13436500	Standard Concentric Glass Nebulizer, K-type, low-Ar-flow, with <i>NebLink</i> sample tube	390
2050-0820N	842-3120-51431	ICAP-63/6500 Concentric Glass Nebulizer low-Ar-flow, with <i>NebLink</i> sample tube	390
TQ-30-A3	TQP-50-A0.5	Concentric Quartz Nebulizer, Type A	590
TR-50-C2	842-3120-51331	<i>Meinhard</i> Type C Concentric Glass Nebulizer, for high salt samples	570
TR-30-C1	472022	'Meinhard' Concentric glass nebulizer, Type C, high solids	390
2050-0920	12698900	TJA Fixed Cross Flow Nebulizer	770
2550-2005	ES-2005-0002	HF-resistant MicroFlow PFA Concentric Nebulizer, 360 ul/min uptake version	990
2550-2028	ME4040	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	570
2550-2024	ME2040	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	990
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
3150-0710	TMM3500	Teflon MiraMist Parallel-Path Nebulizer	990
3150-0720	PMM4000	PEEK MiraMist Parallel-Path Nebulizer	990
2051-0100N		V-Groove High-Solids Glass Nebulizer, with <i>NebLink</i>	780
3150-0730	T2100	Burgener Teflon T2100 'Work-horse' Parallel-Path Nebulizer	990
2028-0340	842-3120-51151	Ball-joint Quartz Adaptor for ICAP 6300/6500 radial	70
2025-0190	13008900	TJA Radial Spray Chamber, Curved Drain, Open Ends	390
2025-0200		TJA Radial Jacketed Spray Chamber, Curved Drain, Open Ends	390
2051-0250	842-3120-51411	Althea Cyclonic Glass Spray Chamber, for ICAP 6300/6500 radial	390
2051-0260	842-3120-51311	Libra Cyclonic Spray Chamber, Baffled, for ICAP 6300/6500 radial	430
2025-0300	395200	TJA Radial, Teflon Baffle for Spray Chamber	120
2052-0300		Teflon Spray Chamber End Cap, for Open-End Spray Chambers, 6mm nebs	215
2027-5180		TJA Radial, MLN Cyclonic Glass Spray Chamber, 12mm drain	450
2027-518J		TJA Radial, Jacketed MLN Cyclonic Glass Spray Chamber, 12mm drain	490
2027-5190	13610500	TJA Radial Cyclonic Spray Chamber, 12mm Drain, 9.6mm Outlet	450
2027-5192	13966300	TJA Radial 'Althea' Cyclonic Spray Chamber with 4mm Drain, 9.6mm Outlet	490
2027-5193		TJA Radial 'Libra' Baffled Cyclonic Spray Chamber, 4mm Drain, 9.6mm Outlet	510
2027-5194		'Electron' Mini-Cyclonic Spray Chamber with 4mm Drain, 9.6mm Outlet	490
2027-5195		'Positron' Baffled Mini-Cyclonic Spray Chamber, 4mm Drain, 9.6mm Outlet	510
2027-519J		TJA Radial Cyclonic Spray Chamber, Water Cooled, 12mm drain	590
EP-3171-3001		PFA Cyclonic spray chamber, O-ring free	990
EP-4157-1000	ES-4157-1000	PC3 Peltier-cooled ESI Quartz Spray Chamber, with Power Supply	2985
EP-9100-5400	ES-9100-5400	PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3970
2025-0400	125407-01	TJA Radial, Quartz 1.0 Injector for Organics	180
2025-0410	125407-00	TJA Radial, Quartz 1.5 Injector for Aqueous	180
2025-0420	136502-01	TJA Radial, Quartz 2.0 Injector	180
2028-0400	842-3120-51341	ICAP-63/6500 Quartz Injector, for organics, 2 Red Rings, ID = 2mm	180
2028-0410	842-3120-51351	ICAP-63/6500 Quartz Injector, Red Ring, ID = 1,5 mm	180
2028-0420	842-3120-51371	ICAP-63/6500 radial Quartz Injector, Blue Ring, ID = 2 mm	180
2028-5090	842-3120-51141	ICAP-6500 radial Quartz Torch	355
2028-5110	842-3120-51741	ICAP-6500 radial Quartz Torch (EMT)	335
2025-0430	135602-02	TJA Radial, Quartz 2.4 Injector	180
2026-0830	136436-01	TJA Radial Torch Alumina Injector, ID = 1,5 mm	570
2026-0840		TJA Radial Torch Alumina Injector, ID = 1,0 mm	570
2027-5060	126432-02	TJA Radial, High/low Flow Torch, 1.5 Aqueous Injector	730
2027-5070	126432-01	TJA Radial, High flow Torch, 1.5 Aqueous Injector	720
2027-507SI	126440-01	TJA Radial, High Flow Torch without Injector	600
2027-507HF	127179-01	TJA Radial, High Flow Torch, 1.5 Alumina Injector	990
2027-5080	125124-01	TJA Radial, High Flow Torch, 1.0 Organics Injector	730
2027-508SI		TJA Radial, High Flow Torch without Injector	600
EP-1423-0100		Quartz 'Drop-In' Injector with ID = 1,0 mm	180
EP-1423-0150		Quartz 'Drop-In' Injector with ID = 1,5 mm	180
EP-1403-0100		Sapphire 'Drop-In' Injector with ID = 1,0 mm	790

EP-1403-0180		Sapphire 'Drop-In' Injector with ID = 1,8 mm	790
EP-1413-0180		Platinum 'Drop-In' Injector with ID = 1,8 mm	1890
EP-1413-0200		Platinum 'Drop-In' Injector with ID = 2,0 mm	1890
2027-5250	085207-00	TJA Torch Anti Arc Cap	55
2027-5260	085207-01	TJA Torch Anti Arc Cap / Flat Edge	60
2027-5350		Argon Connectors for TJA Torches, pkt of 2	30
EPOND No.		Accessory Description	CHF
APEX-Q-TJA		APEX-Q Desolvator, with PolyPro-ST nebulizer	POA
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply	4950
2045-0700	PRG-01	'PERGO' Argon Humidifier, new technology, uses vapour permeation	990
8100-2051	PRL-02	'PEARL' Cleaning Device for concentric glass nebulizers	190
EPME-PG41-100		100-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	360
EPME-PG41-500		500-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	805
		- reference solutions made to measure -	
EPOND N°	OEM P/No.	Description, tubes	CHF
<i>See web site</i>		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12	40
<i>See web site</i>		<i>Cone-terminated</i> Peri-Pump tubes, for low-uptake apps., pkt of 12	48
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40	165
Autosamplers and ICP / AAS Standard and Reference Solutions			
SC-8000-XXXX		Complete range of ESI SC-series autosamplers - please ask for details	POA
ICP-standards		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards	
		- reference solutions made to measure -	
www. e-pond. biz			