

EPOND Part-No.	OEM Part-No.	Description
Nebulizers/Spray chambers/accessories		
Quartz nebulizers		
EP 2050-SILQ-50-05		MEINHARD® SilQ+ Ultra high purity quartz nebulizer
EP 2050-TQP30-C1		MEINHARD® Type C+ concentric quartz nebulizer with PFA
EP 2050-TQP50-A05		MEINHARD® TQ+ quartz nebulizer
EP 2050-MQ30-A2		MEINHARD® Type A quartz nebulizer
Glass nebulizers		
EP 2050-TRP050-05		MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-C2	2010076500	MEINHARD® Type C glass nebulizer, 2ml/min uptake
EP 2050-M030-K2	2010069500	MEINHARD® Type K glass nebulizer, 2ml/min uptake
EP 2050-M030-K3	2010081600	MEINHARD® Type K glass nebulizer, 3ml/min uptake
PFA nebulizers		
ES-2040-28		MicroFlow PFA-ST self-aspirating nebulizer
Poly-Pro nebulizers		
EP 2550-4040-57		Poly-Pro-ST self-aspirating nebulizer
PEEK nebulizer		
EP 3150-0720	G-3161-80000	Burgener PEEK MiraMist nebulizer
Nebulizer cleaner		
PRL02		Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
Nebulizer Gas Quick connects		
EP-F4		Gas Quick connect for MEINHARD® nebulizers with PFA-Press Fit for 4 mm
EP-F4-66		ICP Emission nebulizer Gas Quick Connect with PFA-Press Fit for 4mm for agilent
FG Gas Quick Connect for MEINHARD® nebulizers		
EP-FG		FG Gas Quick Connect for MEINHARD® nebulizers
pergo Argon Nebulizer Gas Humidifier		
PRG-01		pergo Argon high solids without nebulizer, includes gas line connector kit
PRG-01-MaPlus		pergo Argon high solids including MEINHARD® Type A Plus high performance nebulizer
Spray Chambers		
EP 2051-0262	7910043700	Pyrex cyclonic spray chamber for Agilent/Varian 7x0, baffled, 8mm top, Unifit
EP 2080-0190ES		Pyrex cyclonic spray chamber for Agilent/Varian 7x0, baffled, Easy-Seal nebulizer inlet
EP 2080-0210	2010081700	Pyrex cyclonic spray chamber for Agilent/Varian 7x0, no baffle, 8mm top, o-ring inlet
EP 2080-0210ES		Pyrex cyclonic spray chamber for Agilent/Varian 7x0, no baffle, Easy-Seal nebulizer inlet
EP 2080-0220		Pyrex cyclonic spray chamber for Agilent/Varian 7x0, with baffle, o-ring inlet
EP 2080-0130		Pyrex cyclonic spray chamber for Agilent/Varian 7x0, with baffle, Unifit
EP 2080-0170	9910070400	Glass spray chamber, jacketed for coolant
EP 2128-5471		PFA spray chamber with PFA PureCap encap, HF resistant, 47mm
EP 2999-6501	ES-2128-5471	Mounting bracket for 47mm PFA spray chamber
PFA Connector		
ES-4399-2300		Male barbed PFA connector for the transfer line

Torches/shields**Quartz One-Piece Torches**

EP 2080-0120	2010030400	Quartz torch, axial with 90° bend, fixed 2.3mm injector
EP 2080-0500	2010094800	Quartz torch, axial, high solids, fixed injector
EP 2080-0140	2010081200	Quartz torch, axial, fixed 2mm injector
EP 2080-0050	2010069690	Quartz torch, radial, low-flow, fixed 1.5mm injector
EP 2080-0050B	2010069610	Quartz torch with ball joint, radial, fixed 1.5mm injector
EP 2080-0280	2010117400	Quartz torch, radial, fixed 0.8mm injector
EP 2080-0290	2010117500	Quartz torch, radial, fixed 1.4mm injector
EP 2080-0260	2010104700	Quartz torch, axial with right angle end, lowflow 0.8mm I.D. injector
EP 2080-0520	2010095000	Quartz torch, radial, high solids, 1.8mm I.D. injector

Quartz Torches, Demountable Injector

EP 2080-0150	2010083400	Quartz torch, demountable, axial, body only
EP 2080-0100	2010072900	Quartz torch, demountable, radial, body only
EP 2080-1842		Quartz torch, O-ring-free, demountable, use with ESI o-ring free injectors
EP 2080-0480	9910056800	Demountable radial torch kit

Torch with Demountable injector

EP 1842-1100		Quartz torch, demountable, with o-ring-free sapphire injector in PFA base
--------------	--	---

Torch Bonnet

EP 2080-0250	2010070700	Quartz torch bonnet, radial
--------------	------------	-----------------------------

Injectors and accessories**Quartz injectors**

EP 2080-0470	2010077500	Quartz injector, 0.8mm, for organic samples
EP 2080-0460	9910057000	Quartz injector, 1.4mm

Alumina injector

EP 2080-6759	2410067590	Alumina injector, 1.9mm
--------------	------------	-------------------------

O-ring-free Sapphire injector

EP 1865-0180		O-ring-free sapphire injector, 1.8mm
--------------	--	--------------------------------------

O-ring-free Platinum injector

EP 1841-0200		O-ring-free platinum injector, 2mm
--------------	--	------------------------------------

Injector holder

EP 2080-0480	9910057300	PTFE injector holder(Varian-ES series, Liberty, Vista&Ultramass-7000)
--------------	------------	---

Accessories**Copper Load Coils**

EP 2080-0900	110601290	Copper load coil for Agilent/Varian 700-ES, axial
EP 2080-A900	110590490	Copper load coil for Varian Series 2, radial

Unifit adapter

EP 2099-0240	N0776006	Unifit adapter, 6mm glass concentric nebulizer, 12mm inlet port
--------------	----------	---

F0 Encapsulated Micro Tubing

EP-F0-25		F0 Encapsulated micro tubing, 0.25mm I.D., 0.5ml/min, 70cm long, green marker
----------	--	---

F1 Detachable Capillary

EP-F1-80		F1 Detachable Capillary, 0.8mm I.D., 70cm long, blue marker
----------	--	---

F2 Standard PFA Integrated Capillary

EP-F2-50		F2 Standard PFA Integrated Capillary, self-aspiration, 0.5mm I.D., 70cm, orange marker
----------	--	--