

EPOND Part-No.	OEM Part-No.	Description
Nebulizers/Spray chambers/accessories		
Quartz Nebulizers		
EP 2050-SILQ-50-05		MEINHARD® SilQ+Ultra high purity quartz nebulizer
EP 2050-TRP50-A05		MEINHARD® Type TQ+Quartz nebulizer, nominal uptake 0.5 ml
EP 2050-MQ30-A05		MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 30 PSI
EP 2050-MQ30-A1		MEINHARD® Type A quartz nebulizer, 1 ml/min uptake 30 PSI
EP 2050-MQ30-A3		MEINHARD® Type A concentric quartz nebulizer, 3 ml/min uptake
Glass Nebulizers		
EP 2050-TRP50-05		MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A05		MEINHARD® Type A glass nebulizer, 0.5 ml/min uptake
EP 2050-M030-A1	1201318	MEINHARD® Type A glass nebulizer, 1 ml/min uptake, 30 PSI
EP 2050-M030-A3	1040894	MEINHARD® Type A glass nebulizer, 3ml/min uptake
EP 2050-M030-C1		MEINHARD® Type C glass nebulizer, 1ml/min uptake
EP 2050-M030-K1	1040894	MEINHARD® Type K glass nebulizer, 1 ml/min uptake
EP 2050-M030-K2		MEINHARD® Type K glass nebulizer, 2ml/min uptake
PFA Nebulizer		
ES-2040-28		MicroFlow PFA-ST self-aspirating nebulizer
Poly-Pro Nebulizer		
EP 2550-4040-28		Poly-Pro self-aspirating nebulizer
PEEK Nebulizer		
EP 3150-0720		PEEK MiraMist nebulizer
pergo Argon Nebulizer Gas Humidifier		
PRG-01		<i>pergo</i> Argon high solids without nebulizer, includes gas line connector kit
PRG-01-MaPlus		<i>pergo</i> Argon high solids including MEINHARD® Type A Plus high performance nebulizer
Nebulizer Cleaner		
PRL02	8100-2051	Pearl Nebulizer Cleaner for all Quartz and glass concentric nebulizers
Conical Spray Chamber		
EP 2045-0220	3204313	Pyrex conical spray chamber, impact, jacketed
EP 2045-0230	3600374	Quartz conical spray chamber, impact, jacketed
EP 2045-0810	3600170	Quartz conical spray chamber, Peltier cooled
Scott Spray Chamber		
EP 2045-0160	3200841	Quartz Scott spray chamber, jacketed for coolant
PFA Spray Chamber		
EP 2124-0351		35mm PFA spry chamber and PFA Purecap Endcap for PQ2, PQ3&Axiom
Spray Chamber Cooling Block&Mounting Brackets		
EP 6124-0001		Spray chamber cooling block with mounting bracket, anodized aluminum
EP 6124-1002		Mounting bracket for aluminum spray chamber cooling block for VG PQ2
EP 6124-2002		Mounting bracket for aluminum spray chamber cooling block for VG Axiom
Torches/shields		
Torch Shields/Screens		
EP 2045-9001	3600157	Platinum torch shield
EP 2045-9900	3601219	Silver torch screen

		Quartz One-Piece Torches with Ball Joint
EP 2045-A050	3201192	Quartz torch with ball joint, axial, fixed 1.5mm injector
EP 2045-0050	319-308	Quartz torch with o-ring seal ball joint, axial, fixed 1.5mm injector
EP 2045-0060	3201732	Quartz torch with ball joint, axial, fixed 2mm injector
		Quartz Torches, Demountable Injector
EP 2045-0080	1200080	Quartz torch, demountable, axial
EP 1202-0201		ESI precision o-ring-free quartz torch, demountable, axial
		Torch Bonnet
EP 2045-0105	319-311	Quartz torch bonnet for Thermo PlasmaQuad/NU
		Torch Shield
EP 2045-0290	3203789	Quartz torch shield for Thermo PlasmaQuad
	3600884	

Cones

		Sampling Cones
VG1001-Ni	3004661	Ni sample cone
VG1021-Ni	3600812	Ni sample cone, Nicone TM, Cu base
EP 2045-0485	1600485	Ni sample cone(HPI) supercone
VG1026-Pt	1600335	Pt sample cone(HP) supercone
EP 2045-0521	3202521	Pt sampler cone(6.35mm) HPI
EP 2045-0500	1067500	Pt sampler cone, Ni core
ES-3000-1808	1067501	Pt sampler, Ni core
EP-01289	3601289	Pt sampler cone
Nu1001-GG	3004382	Gasket, graphite for sampler, pack of 3
EP 2048-9000	1310900	Gasket, graphite for sampler, holes, 3 pk
		Skimmer Cones
VG1004-Ni	3200860	Ni micro-skimmer cone
VG1022-Ni	3600811	Ni skimmer cone, Nicone TM
VG1008-Pt	3201101	Pt micro-skimmer cone, XS
EP-47460	1047460	Pt skimmer, Ni base(6mm)HPI
VG1023-00	3600813	Adapter, copper skimmer with screws
EP-01240	120124	O-ring, Viton for micro-skimmer, pack of 3
		Cone Cleaning Kit&Detergent
EPMNTKIT1		Cone Cleaning Kit
EPCitranox		Citranox Detergent

Accessories

		Quartz Elbows
EP 2045-0300	3600169	Quartz elbow, ball/socket 17°
EP 2045-0320	3200727	Quartz elbow, ball/socket 90° for Scott type quartz spray chamber
		Copper Load Coil
EP 2045-0618		Copper load coil for Thermo PQ2, PQ3, Excell
		Electron Multipliers
EPMAF214	2203433	ETP Electron Multiplier for Thermo PQ3 Excel(simultaneous)
EPMAF562A	2203405	ETP Electron Multiplier for VG PQ3 Sequential(SXP Rods)