

Precision Quartz Torches and O-ring Free Injectors

Part of a complete o-ring free sample introduction system

- ▲ Direct replacement for standard torches and injectors
- ▲ O-ring-free Teflon base
 - reduces contamination
 - easy to clean
 - easy to insert
- ▲ Precision alignment of torch and injector eliminates injector wobble and off-axis deflection.
- ▲ Precision o-ring-free injector with PFA base available in a variety of sizes and materials including:
 - Quartz
 - Sapphire
 - Platinum
- ▲ Stronger, reinforced torch base helps prevent accidental breakage during maintenance
- ▲ Available for most ICP and ICPMS instruments



The precision quartz torch and o-ring free injector directly replace the standard demountable torch and injector.